

---

**Anti-EphA1**

Cat #: HM1128  
Rabbit polyclonal IgG  
0.2 µg/µl, store at 4 °C

For research use only

**BACKGROUND**

The ephrin receptors are the largest group of receptor protein tyrosine kinases. They are implicated in mediating numerous developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands.

**SPECIFICITY**

This antibody specifically recognizes EphA1 of human, mouse and rat origin and can be used in Western blotting, immunoprecipitation and immunostaining.

**IMMUNOGEN**

A peptide mapping at the carboxy terminus of human EphA1.

**STORAGE**

This antibody is stable for 12 months when stored at 2-8°C.

**REFERENCES**

1. Beckmann, M.P., Cerretti, D.P., Baum, P., VandenBos, T., James, L., Farrah, T., Kozlosky, C., Hollingsworth, T., Shilling, H., Maraskovsky, E., Fletcher, F.A., Lhotak, V., Pawson, T., and Lyman, S.D. 1994. Molecular characterization of a family of ligands for eph-related tyrosine kinase receptors. *EMBO J.* 13: 3757-3762.
2. Cheng, H.-J. and Flanagan, J.G. 1994. Identification and cloning of ELF-1, a developmentally expressed ligand for the Mek4 and Sek receptor tyrosine kinases. *Cell* 79: 157-168.
3. Fox, G.M., Holst, P.L., Chute, H.T., Lindberg, R.A., Janssen, A.M., Basu, R., and Welcher, A.A. 1995. DNA cloning and tissue distribution of five human EPH-like receptor proteintyrosine kinases. *Oncogene* 10: 897-905.
4. Ciossek, T., Millauer, B., and Ullrich, A. 1995. Identification of alternatively spliced mRNAs encoding variants of MDK1, a novel receptor tyrosine kinase expressed in the murine nervous system. *Oncogene* 10: 97-108.
5. Valenzuela, D.M., Rojas, E., Griffiths, J.A., Compton, D.L., Gisser, M., Ip, N.Y., Goldfarb, M., and Yancopoulos, G.D. 1995. Identification of full length and truncated forms of Ehk-3, a novel member of the Eph receptor tyrosine kinase family. *Oncogene* 10: 1573-1580.
6. Kozlosky, C.J., Eugene Maraskovsky, E., McGrew, J.T., VandenBos, T., Teepe, M., Lyman, S.D., Srinivasan, S., Fletcher, F.A., Gayle, R.B. III., Cerretti, D.P., and Beckmann, M.P. 1995. Ligands for the receptor tyrosine kinases hek and elk: isolation of cDNAs encoding a family of proteins. *Oncogene* 10: 299-306.

7. Bennett, B.D., Zeigler, F.C., Gu, Q., Fendly, B., Goddard, A.D., Gillett, N., and Matthews, W. 1995. Molecular cloning of a ligand for the EPH-related receptor protein-tyrosine kinase Htk. *Proc. Natl. Acad. Sci. USA* 92: 1866-1870.
8. Gale, N.W., Holland, S.J., Valenzuela, D.M., Flenniken, A., Pan, L., Ryan, T.E., Henkemeyer, M., Strebhardt, K., Hirai, H., Wilkinson, D.G., Pawson, T., Davis, S. and Yancopoulos, G.D. (1996) Eph receptors and ligands comprise two major specificity subclasses and are reciprocally compartmentalized during embryogenesis. *Neuron* 17, 9-19.
9. Wilkinson, D.G. (2001) Multiple roles of EPH receptors and ephrins in neural development. *Mol. Cell. Probes* 2, 155-164.

**PRODUCT FROM HYPROMATRIX, INC.****A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000

**B. Staining AntibodyArray™s**

1. Staining AntibodyArray™ I  
Catalog Number HM8100
2. AntibodyArray Staining Apparatus  
Catalog Number HM8000

**C. Antibodies****1. HRP-conjugated antibodies**

- anti-phosphotyrosine  
Catalog Number HM2040
- anti-phosphoserine  
Catalog Number HM2070
- anti-phosphothreonine  
Catalog Number HM2090

and more...

**2. Primary antibodies**

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at [www.hypromatrix.com](http://www.hypromatrix.com).

**CONTACT**

**Hypromatrix, Inc.**  
100 Barber Avenue  
Worcester, MA 01606  
USA

Tel: 508-856-7900  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)  
Web: [www.hypromatrix.com](http://www.hypromatrix.com)