
anti-BMX

Cat #: HM1036
Goat polyclonal IgG
0.2 µg/µl, store at 4 °C

For research use only

BACKGROUND

BMX is a member of the Btk family of non-receptor tyrosine kinases. Other members of this family include Tec, Emt (also known as Itk or Tsk), Btk (previously known as Atk), BPK(Emb), Txk (also known as Rlk) and Dsrc28C. Most members of the family contain SH3 and SH2 domains, a pleckstrin homology (PH) and a Tec homology (TH) domain. The 80 kDa Bmx protein is expressed at highest levels in the heart.

SPECIFICITY

This antibody specifically recognizes BMX of human, mouse and rat origin.

The antibody can be used in Western blotting, immunoprecipitation and immunostaining.

IMMUNOGEN

A recombinant peptide corresponding to the C terminus of human BMX.

STORAGE

This antibody is stable for 12 months when stored at 2-8°C.

REFERENCE

1. Tamagnone, L., Lahtinen, I., Mustonen, T., Virtaneva, K., Francis, F., Muscatelli, F., Alitalo, R., Smith, C.I.E., Larsson, C. and Alitalo, K. 1994. BMX, a novel nonreceptor tyrosine kinase gene of the BTK/ITK/TEC/TXK family located in chromosome Xp22.2. *Oncogene* 9: 3683-3688.
2. Yamada, N., Kawakami, Y., Kimura, H., Fukamachi, H., Baier, G., Altman, A., Kato, T., Inagaki, Y., and Kawakami, T. 1993. Structure and expression of novel protein tyrosine kinases, Emb and Emt, in hematopoietic cells. *Biochem. Biophys. Res. Comm.* 192: 231-240.
3. Haire, R.N., Ohta, Y., Lewis, J.E., Fu, S.M., Kroisel, P., and Litman, G.W. 1994. TXK, a novel human tyrosine kinase expressed in T cells shares sequence identity with Tec family kinases and maps to 4p12. *Hum. Mol. Genet.* 3: 897-901.
4. Saharinen P, Ekman N, Sarvas K, Parker P, Alitalo K, Silvennoinen O. (1997) The Bmx tyrosine kinase induces activation of the Stat signaling pathway, which is specifically inhibited by protein kinase Cdelta. *Blood.* 90:4341-53.
5. Qiu Y, Robinson D, Pretlow TG, Kung HJ. (1998) Etk/Bmx, a tyrosine kinase with a pleckstrin-homology domain, is an effector of phosphatidylinositol 3'-kinase and is involved in interleukin 6-induced neuroendocrine differentiation of prostate cancer cells. *Proc Natl Acad Sci U S A.* 95:3644-9.
6. Mao J, Xie W, Yuan H, Simon MI, Mano H, Wu D. (1998) Tec/Bmx non-receptor tyrosine kinases are involved in

regulation of Rho and serum response factor by Galpha12/13. *EMBO J.* 17:5638-46.

7. Ekman N, Arighi E, Rajantie I, Saharinen P, Ristimaki A, Silvennoinen O, Alitalo K. (2000) The Bmx tyrosine kinase is activated by IL-3 and G-CSF in a PI-3K dependent manner. *Oncogene.* 19:4151-8.
8. Wu YM, Huang CL, Kung HJ, Huang CY. (2001) Proteolytic activation of ETK/Bmx tyrosine kinase by caspases. *J Biol Chem.* 276:17672-8.
9. Bagheri-Yarmand R, Mandal M, Taludker AH, Wang RA, Vadlamudi RK, Kung HJ, Kumar R. (2001) Etk/Bmx tyrosine kinase activates Pak1 and regulates tumorigenicity of breast cancer cells. *J Biol Chem.* 276:29403-9.

PRODUCTS FROM HYPROMATRIX, INC.**A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™
Catalog Number HM3000
2. Apoptosis AntibodyArray™
Catalog Number HM4000
3. Cell Cycle AntibodyArray™
Catalog Number HM5000

B. Staining AntibodyArray™s

1. Staining AntibodyArray™ I
Catalog Number HM8100
2. AntibodyArray Staining Apparatus
Catalog Number HM8000

C. Antibodies**1. HRP-conjugated antibodies**

- anti-phosphotyrosine
Catalog Number HM2040
- anti-phosphoserine
Catalog Number HM2070
- anti-phosphothreonine
Catalog Number HM2090

and more...

2. Primary antibodies

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at www.hypromatrix.com.

CONTACT

Hypromatrix, Inc.
100 Barber Avenue
Worcester, MA 01606
USA

Tel: 508-856-7900
Fax: 508-302-0748
Email: contact@hypromatrix.com
Web: www.hypromatrix.com