

---

**anti-c-IAP1**

Cat #: HM1087  
Rabbit polyclonal IgG  
0.2 µg/µl, store at 4 °C

For research use only

**BACKGROUND**

C-IAP1 is a mammalian homologue of baculovirus protein p35. It is a member of a family of proteins that inhibits apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. C-IAP1 inhibits apoptosis induced by serum deprivation but does not affect apoptosis resulting from exposure to menadione, a potent inducer of free radicals.

**SPECIFICITY**

This antibody specifically recognizes c-IAP1 of human, mouse and rat origins.

The antibody can be used in Western blotting, immunoprecipitation and immunostaining.

**IMMUNOGEN**

Full-length recombinant human C-IAP1 protein.

**REFERENCES**

1. Rothe,M., Pan,M.G., Henzel,W.J., Ayres,T.M. and Goeddel,D.V. (1995) The TNFR2-TRAF signaling complex contains two novel proteins related to baculoviral inhibitor of apoptosis proteins. *Cell* 83, 1243-1252.
2. Liston,P., Roy,N., Tamai,K., Lefebvre,C., Baird,S., Cherton-Horvat,G., Farahani,R., McLean,M., Ikeda,J.E., MacKenzie,A. and Korneluk,R.G. (1996) Suppression of apoptosis in mammalian cells by NAIP and a related family of IAP genes. *Nature* 379, 349-353.
3. Deveraux QL, Takahashi R, Salvesen GS, Reed JC. X-linked IAP is a direct inhibitor of cell-death proteases. *Nature*. 1997 Jul 17;388:300-4.
4. Verhagen AM, Ekert PG, Pakusch M, Silke J, Connolly LM, Reid GE, Moritz RL, Simpson RJ, Vaux DL. Identification of DIABLO, a mammalian protein that promotes apoptosis by binding to and antagonizing IAP proteins. *Cell*. 2000 Jul 7;102:43-53.
5. Wang,Q., Wang,X. and Evers,B.M. (2003) Induction of cIAP-2 in human colon cancer cells through PKC delta/NF-kappa B. *J. Biol. Chem.* 278, 51091-51099.
6. Yang,Q.H. and Du,C. (2004) Smac/DIABLO selectively reduces the levels of c-IAP1 and c-IAP2 but not that of XIAP and livin in HeLa cells. *J. Biol. Chem.* 279, 16963-16970.
7. Hu,S. and Yang,X. (2003) Cellular inhibitor of apoptosis 1 and 2 are ubiquitin ligases for the apoptosis inducer Smac/DIABLO. *J. Biol. Chem.* 278, 10055-10060.
8. Gagnon,V., St-Germain,M.E., Parent,S. and Asselin,E. (2003) Akt activity in endometrial cancer cells: regulation of cell survival through cIAP-1. *Int. J. Oncol.* 23, 803-810.
9. Plenchette,S., Cathelin,S., Rebe,C., Launay,S., Ladoire,S., Sordet,O., Ponnelle,T., Debili,N., Phan,T.H., Padua,R.A., Dubrez-Daloz,L. and Solary,E. (2004) Translocation of the inhibitor of apoptosis protein c-IAP1 from the nucleus to the Golgi in hematopoietic cells undergoing differentiation: a nuclear export signal-mediated event. *Blood* 104 (7), 2035-2043

**PRODUCT FROM HYPROMATRIX, INC.****A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000

**B. Staining AntibodyArray™s**

1. Staining AntibodyArray™ I  
Catalog Number HM8100
2. AntibodyArray Staining Apparatus  
Catalog Number HM8000

**C. Antibodies****1. HRP-conjugated antibodies**

- anti-phosphotyrosine  
Catalog Number HM2040
- anti-phosphoserine  
Catalog Number HM2070
- anti-phosphothreonine  
Catalog Number HM2090

and more...

**2. Primary antibodies**

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at [www.hypromatrix.com](http://www.hypromatrix.com).

**CONTACT**

**Hypromatrix, Inc.**  
100 Barber Avenue  
Worcester, MA 01606  
USA

Tel: 508-856-7900  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)  
Web: [www.hypromatrix.com](http://www.hypromatrix.com)