

---

**Anti-TR**

Cat #: HM1370  
Rabbit polyclonal IgG  
0.2 µg/µl, store at 4 °C

For research use only

**BACKGROUND**

Thyroid hormone nuclear receptors (TRs) are ligand-dependent transcription factors which regulate growth, differentiation and development. They have been shown to mediate the biological activities of thyroid hormone. TRs bind to thyroid hormone response elements (TREs) with half site binding motifs in the orientation of palindrome, direct repeat or inverted palindrome. The genes encoding TR α1 and TR β1 are mapped to human chromosomes 17 and 3, respectively. Knockout studies in mice suggest that the different receptors, while having certain extent of redundancy, may mediate different functions of thyroid hormone.

**SPECIFICITY**

This antibody specifically reacts with TR α1 and β 1 of human, mouse and rat origin.

The antibody can be used in Western blotting, immunoprecipitation and immunohistochemistry.

**IMMUNOGEN**

A recombinant protein corresponding to the full-length human thyroid hormone receptor α 1.

**STORAGE**

This antibody is stable for 12 months when stored at 2-8° C.

**REFERENCES**

1. Mangelsdorf, D.J., Umesono, K., and Evans, R.M. 1994. In Sporn, M.B., Roberts, A.B., and Goodman, D.S., eds. THE RETINOIDS: Biology, Chemistry, and Medicine, 2nd Edition. New York: Raven Press, Ltd.: 319-349.
2. Näär, A., Boutin, J.M., Lipkin, S.M., Yu, V.C., Holloway, J.M., Glass, C.K., and Rosenfeld, M.G. 1991. The orientation and spacing of core DNA-binding motifs dictate selective transcriptional responses to three nuclear receptors. *Cell* 65: 1267-1271.
3. Lazar, M.A. 1993. Thyroid hormone receptors: multiple forms, multiple possibilities. *Endocrinol. Rev.* 14: 184-193.
4. Meier, C.A., Parkison, C., Chen, A., Ashizawa, K., Meier-Heusler, S.C., Muchmore, P., Cheng, S.Y., and Weintraub, B.D. 1993. Interaction of human beta 1 thyroid hormone receptor and its mutants with DNA and retinoid X receptor beta. T3 response element-dependent dominant negative potency. *J. Clin. Invest.* 92: 1986-1993.
5. Zhang, X.K. and Pfahl, M. 1993. Hetero- and homodimeric receptors in thyroid hormone and vitamin A action. *Receptor* 3: 183-191.
6. Bhat, M.K., Ashizawa, K., and Cheng, S.Y. 1994. Phosphorylation enhances the target gene sequence-dependent dimerization of thyroid hormone receptor with

retinoid X receptor. *Proc. Natl. Acad. Sci. USA* 91: 7927-7931.

7. Sugawara, A., Yen, P.M., Apriletti, J.W., Ribeiro, R.C.J., Sacks, D.B., Baxter, J.D., and Chin, W.W. 1994. Phosphorylation selectively increases triiodothyronine receptor homodimer binding to DNA. *J. Biol. Chem.* 269: 433-437.
8. Forrest, D., Reh, T.A. and Rusch, A. (2002) Neurodevelopmental control by thyroid hormone receptors. *Curr. Opin. Neurobiol.* 12, 49-56.
9. Lin, H.M., Zhao, L. and Cheng, S.Y. (2002) Cyclin D1 Is a Ligand-independent Co-repressor for Thyroid Hormone Receptors. *J. Biol. Chem.* 277, 28733-28741.

**PRODUCTS FROM HYPROMATRIX, INC.****A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000

**B. Staining AntibodyArray™s**

1. Staining AntibodyArray™ I  
Catalog Number HM8100
2. AntibodyArray Staining Apparatus  
Catalog Number HM8000

**C. Antibodies****1. HRP-conjugated antibodies**

- anti-phosphotyrosine  
Catalog Number HM2040
- anti-phosphoserine  
Catalog Number HM2070
- anti-phosphothreonine  
Catalog Number HM2090

and more...

**2. Primary antibodies**

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at [www.hypromatrix.com](http://www.hypromatrix.com).

**CONTACT**

**Hypromatrix, Inc.**  
100 Barber Avenue  
Worcester, MA 01606  
USA

Tel: 508-856-7900  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)  
Web: [www.hypromatrix.com](http://www.hypromatrix.com)