

---

**anti-VDAC1**

Cat #: HM1391  
Goat polyclonal IgG  
0.2 µg/µl, store at 4 °C

For research use only

**BACKGROUND**

The Voltage-Dependent Anion Channels (VDAC) are outer membrane mitochondrial proteins. They form the major pores through which adenine nucleotides cross the outer mitochondrial membrane. Together with adenine nucleotide translocator (ANT) VDAC1 and VDAC2 form the permeability transition pore complex (PTPC) of the mitochondria. The complex is thought to allow the mitochondria to undergo metabolic uncoupling and irreversible morphologic changes that ultimately destroy the mitochondria during apoptosis.

**SPECIFICITY**

This antibody specifically reacts with VDAC1 of human, mouse and rat origin.

The antibody can be used in Western blotting, immunoprecipitation and immunohistochemistry.

**IMMUNOGEN**

A synthetic peptide derived from N-terminus of human VDAC1 protein.

**STORAGE**

This antibody is stable for 12 months when stored at 2-8°C.

**REFERENCES**

1. Cozens, A.L., Runswick, M.J., and Walker, J.E. 1989. DNA sequences of two expressed nuclear genes for human mitochondrial ADP/ATP translocase. *J. Mol. Biol.* 206: 261-280.
2. Li, K., Warner, C.K., Hodge, J.A., Minoshima, S., Kudoh, J., Fukuyama, R., Maekawa, M., Shimizu, Y., Shimizu, N., and Wallace, D.C. 1989. A human muscle adenine nucleotide translocator gene has four exons, is located on chromosome 4, and is differentially expressed. *J. Biol. Chem.* 264: 13998-14004.
3. Blachly-Dyson, E., Zambronicz, E.B., Yu, W.H., Adams, V., McCabe, E.R., Adelman, J., Colombini, M., and Forte, M. 1993. Cloning and functional expression in yeast of two human isoforms of the outer mitochondrial membrane channel, the voltagedependent anion channel. *J. Biol. Chem.* 268: 1835-1841.
4. Green, D.R. and Reed, J.C. 1998. Mitochondria and apoptosis. *Science* 281: 1309-1312.
5. Zamzami, N., Susin, S.A., Marchetti, P., Hirsch, T., Gomez-Monterrey, I., Castedo, M., and Kroemer, G. 1996. Mitochondrial control of nuclear apoptosis. *J. Exp. Med.* 183: 1533-1544.
6. Marzo, I., Brenner, C., Zamzami, N., Jurgens-meier, J.M., Susin, S.A., Viera, H.L.A., Prevost, M.C., Xie, Z., Matsuyama, S., Reed JC, and Kroemer, G. 1998. Bax and adenine nucleotide translocator cooperate in the mitochondrial control of apoptosis. *Science* 281: 2027-2031.
7. Sugiyama, T., Shimizu, S., Matsuoka, Y., Yoneda, Y. and Tsujimoto, Y. (2002) Activation of mitochondrial voltage-dependent anion channel by a pro-apoptotic BH3-only protein Bim. *Oncogene* 21, 4944-4956.

**PRODUCT FROM HYPROMATRIX, INC.****A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000

**B. Staining AntibodyArray™s**

1. Staining AntibodyArray™ I  
Catalog Number HM8100
2. AntibodyArray Staining Apparatus  
Catalog Number HM8000

**C. Antibodies****1. HRP-conjugated antibodies**

- anti-phosphotyrosine  
Catalog Number HM2040
- anti-phosphoserine  
Catalog Number HM2070
- anti-phosphothreonine  
Catalog Number HM2090

and more...

**2. Primary antibodies**

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at [www.hypromatrix.com](http://www.hypromatrix.com).

**CONTACT**

**Hypromatrix, Inc.**  
100 Barber Avenue  
Worcester, MA 01606  
USA

Tel: 508-856-7900  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)  
Web: [www.hypromatrix.com](http://www.hypromatrix.com)