

---

**Anti- Sp2**

Cat #: HM1341  
Rabbit polyclonal IgG  
0.2 µg/µl, store at 4 °C

For research use only

**BACKGROUND**

The Sp zinc finger transcription factor family is composed of four members (Sp1, Sp2, Sp3 and Sp4) that share similar structures. All four proteins contain a highly conserved DNA-binding domain composed of three zinc fingers at the C-terminus. They bind the consensus sequence GGGGCGGGC and other closely related sequences which are known as GC boxes. Sp1, Sp3 and Sp4 are more closely related to each other than to Sp2, which does not bind to a GC-box but to a GT-rich element. Sp1, Sp2 and Sp3 are ubiquitously expressed, while Sp4 is abundantly expressed in brain with limited expression in other tissues.

**SPECIFICITY**

This antibody specifically reacts with Sp2 of human, mouse and rat origin.

The antibody can be used in Western blotting, immunoprecipitation and immunohistochemistry.

Molecular Weight of Sp2: 80 kDa.

**IMMUNOGEN**

A synthetic peptide derived from N-terminus of human Sp2 protein.

**STORAGE**

This antibody is stable for 12 months when stored at 2-8°C.

**REFERENCES**

1. Dynan, W.S. et al. 1983. Isolation of transcription factors that discriminate between different promoters recognized by RNA polymerase II. *Cell* 32: 669-680.
2. Dynan, W.S. et al. 1983. The promoter-specific transcription factor Sp1 binds to upstream sequences in the SV40 early promoter. *Cell* 35: 79-87.
3. Kadonaga, J.T., et al. 1987. Isolation of cDNA encoding transcription factor Sp1 and functional analysis of the DNA binding domain. *Cell* 51: 1079-1090.
4. Kadonaga, J.T., Courey, A.J., Ladika, J., and Tjian, R. 1988. Promoter-selective activation of transcription by Sp1. in *The Control of Human Retrovirus Gene Expression*. Franza, B.R. Jr., Cullen, B.R., and Wong-Staal, F., eds. Cold Spring Harbor Laboratory, Cold Spring Harbor, NY: 239-250.
5. Kingsley, C. and Winoto, A. 1992. Cloning of GT box-binding proteins: a novel Sp1 multigene family regulating T-cell receptor gene expression. *Mol. Cell. Biol.* 12: 4251-4261.
6. Hagen, G., Müller, S., Beato, M., and Suske, G. 1992. Cloning by recognition site screening of two novel GT box

binding proteins: a family of Sp1 related genes. *Nucleic Acids Res.* 20: 5519-5525.

7. Hagen, G., Müller, S., Beato, M., and Suske, G. 1994. Sp1-mediated transcriptional activation is repressed by Sp3. *EMBO J.* 13: 3843-3851.
8. Hagen, G., Dennig, J., Preiss, A., Beato, M., and Suske, G. 1995. Functional analyses of the transcription factor Sp4 reveal properties distinct from Sp1 and Sp3. *J. Biol. Chem.* 270: 24989-24994.
9. Bigger, C.B., Melnikova, I.N., and Gardner, P.D. 1997. Sp1 and Sp3 regulate expression of the neuronal nicotinic acetylcholine receptor beta4 subunit gene. *J. Biol. Chem.* 272: 25976-25982.

**PRODUCTS FROM HYPROMATRIX, INC.****A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000

**B. Staining AntibodyArray™s**

1. Staining AntibodyArray™ I  
Catalog Number HM8100
2. AntibodyArray Staining Apparatus  
Catalog Number HM8000

**C. Antibodies****1. HRP-conjugated antibodies**

- anti-phosphotyrosine  
Catalog Number HM2040
- anti-phosphoserine  
Catalog Number HM2070
- anti-phosphothreonine  
Catalog Number HM2090

and more...

**2. Primary antibodies**

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at [www.hypromatrix.com](http://www.hypromatrix.com).

**CONTACT**

**Hypromatrix, Inc.**  
100 Barber Avenue  
Worcester, MA 01606  
USA

Tel: 508-856-7900  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)

