

## Custom AntibodyArray™

Cat #: HM6000

### APPLICATIONS

Each Custom AntibodyArray™ contains high quality antibodies against proteins of customers' choice. Custom AntibodyArray™ systems are generally used in the following applications:

To confirm or further study the positive results from other types of AntibodyArray™ assays.

### PRODUCT

Each Custom AntibodyArray™ contains antibodies immobilized on a nitrocellulose membrane. The list of the antibodies on the array is included in the Instruction Manual.

### STORAGE

Store at 4°C, do not freeze; stable for 6 months from the date of shipment.

### RESEARCH USE

For research use only, not for use in diagnostic or therapeutic procedures.

### REFERENCE

Wang Y, Wu TR, Cai S, Welte T, Chin YE. 2000. Stat1 as a component of tumor necrosis factor alpha receptor 1-TRADD signaling complex to inhibit NF-kappaB activation. Mol. Cell. Biol. 20(13):4505-4512

## PRODUCT FROM HYPROMATRIX, INC.

### A. AntibodyArray™s:

1. Signal Transduction AntibodyArray™  
Catalog Number HM3000
2. Apoptosis AntibodyArray™  
Catalog Number HM4000
3. Cell Cycle AntibodyArray™  
Catalog Number HM5000
4. Custom AntibodyArray™  
Catalog Number HM6000

### B. Antibodies

#### 1. HRP-conjugated antibodies

- anti-GST  
Catalog Number HM2010
- anti-GFP  
Catalog Number HM2020
- anti-c-Myc  
Catalog Number HM2030
- anti-phosphotyrosine  
Catalog Number HM2040

and more...

#### 2. Primary antibodies

Hypromatrix offers a variety of high quality antibodies to our array customers. For a complete list of antibodies and their specificities, please contact Hypromatrix.

### CONTACT

**Hypromatrix, Inc.**  
25 Winthrop St., C/O MBI  
Worcester, MA 01604  
USA

Tel: 508-797-1700  
Fax: 508-302-0748  
Email: [contact@hypromatrix.com](mailto:contact@hypromatrix.com)  
Web: [www.hypromatrix.com](http://www.hypromatrix.com)